Cromatógrafos SHIMADZU

Dr. Marcelo Del Grande
Sinc do Brasil
mdelgrande@sinc.com.br

GC 2010

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



5. Automatic Adjustment of Retention Time (AART)

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
Automatic Adjustment of Retention Time (A


Removing the contaminated part of capillary column for maintenance.
Column length is short


The retention time becomes shorter.

The data can not be processed properly.

Automatic Adjustment of Retention Time (AART)


SHIMADZU AART function


AART : Automatic Adjustment of Retention time

Principle of Retention Time Adjustment

Method information

After cutting the column

□ □ Target compound

Libraries Options

1. NIST/EPA – 190.825 compounds
2. Wiley – 338.323 compounds
3. Pesticides – 318 EI and 173 CI/NCI compounds
4. Drogas – 4.366 compounds
5. Customized
- 6. FFNSC – 1.156 compounds (Exclusive library)**

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

Aplicações

Petrochemical

- Aroma in gasoline
- Alcohols in gasoline
- Hydrocarbons


Flavor&Fragrances and Natural Products

- Essential oil
- Fast analysis

Drugs

- Benzodiazepines

Analysis of Lime Oil


ID	RT	RI	Compound Name
1	7.367	926	alpha- Thujene
2	7.608	933	alpha- Pinene
3	8.142	947	Camphene
4	9.058	973	Sabinene
5	9.208	977	beta- Pinene
6	9.708	991	Myrcene
7	10.758	1017	alpha- Terpinene
8	11.108	1025	Cymene
9	11.408	1032	Limonene

ID	RT	RI	Compound Name
10	12.6	1059	gamma- Terpinene
11	13.892	1089	Terpinolene
12	18.65	1192	alpha- Terpineol
13	21.025	1243	Citral
15	26.667	1366	Neryl acetate
16	27.525	1385	Geranyl acetate
17	29.167	1423	Valencene
18	29.833	1438	alpha- trans- Bergamotene
19	32.95	1512	beta- Bisabolene

Identification Using FFNSC Library in Fast-GC/MS

**Sample :Lavender
 Essential Oil**

In case of Fast-GC/MS, The retention index shifts approximately 20 comparing with conventional GC/MS.

The retention index allowance in fast-GC/MS is needed to set more extensive than in conventional-GC/MS.

ID	Compound Name	Fast				Conventional
		RT	Actual RI	Register RI	Error	Error
1	alpha- Pinene	1.577	938	933	5	-2
2	Camphene	1.677	954	953	1	-6
3	beta- Pinene	1.855	983	978	5	-3
4	Myrcene	1.922	993	991	2	-1
5	Cymene	2.160	1030	1024	6	0
6	Limonene	2.188	1034	1030	4	-2
7	Eucalyptol	2.217	1039	1032	7	0
8	gamma- Terpinene	2.385	1064	1058	6	0
9	Linalool	2.668	1108	1101	7	3
10	Camphor	2.973	1156	1149	7	-3
11	Isoborneol	3.095	1175	1165	10	-8
12	alpha- Terpeneol	3.240	1199	1198	1	-7
13	Linalyl acetate	3.600	1260	1250	10	7
14	Bornyl acetate	3.800	1294	1285	9	2
15	Neryl acetate	4.212	1368	1361	7	4
16	(Z)- beta- Farnesene	4.703	1462	1454	8	3
17	Germacrene D	4.887	1497	1480	17	4
18	beta- Bisabolene	4.988	1518	1508	10	3
19	Caryophyllene oxide	5.408	1605	1587	18	1


Benzodiazepines and Metabolites


Alprazolam <i>α</i> -Hydroxyalprazolam	Estazolam M-□ M-□	Medazepam Desmethyldiazepam Oxazepam
Bromazepam 3-OH ABBP	Etizolam 1-Hydroxyetizolam <i>α</i> -Hydroxyetizolam	Mexazolam Delorazepam
Brotizolam a-Hydroxybrotizolam 4-Hydroxybrotizolam	Fludiazepam Desalkylflurazepam	Nimetazepam 7-Aminonimetazepam 7-Acetamidonimetazepam
Chlordiazepoxide Desmethyldiazepam Oxazepam	Flutazolam 1-Ethanolflurazepam	Nitrazepam 7-Aminonitrazepam
Clobazam N-Desmethyloclobazam	Flutoprazepam Desalkylflurazepam	Oxazepam Oxazolam Desmethyldiazepam Oxazepam
Clonazepam 7-Aminoclonazepam 7-Acetamidoclonazepam	Flunitrazepam [3-Hydroxyflunitrazepam] 7-Aminoflunitrazepam [7-Acetamidoflunitrazepam]	Prazepam 3-Hydroxyprazepam Oxazepam
Clorazepic acid Desmethyldiazepam Oxazepam	Ethyl loflazepate Desalkylflurazepam	Quazepam Desmethyldiazepam M-4
Clotiazepam Desmethylclotiazepam Y-10247	Flurazepam 1-Ethanolflurazepam Desalkylflurazepam	Temazepam Oxazepam
Delorazepam	[Haloxazolam]	Triazolam <i>α</i> -Hydroxytriazolam
Diazepam Desmethyldiazepam Oxazepam	RAZ-609 Lorazepam	

Black: Benzodiazepine Blue: Benzodiazepine metabolites [] shows the compounds which are not included in this study due to low sensitivity


Special Configuration

Comprehensive GCXGCMSq (MDGCMS)

MDGCMS


Transferred Chiral Compounds

- 1 - α -pinene
- 2 - sabinene
- 3 - β -pinene
- 4 - limonene
- 5 - linalool
- 6 - terpinen-4-ol
- 7 - α -terpineol
- 8 - linalyl acetate


1st detector of bergamot oil

MDGCMS

MDGCMS

2nd detector of bergamot oil

Summary Applications

- It is important to use an appropriate system for getting efficient data when you use GC/MS.

	Environmental	Food	Medical / Pharmaceutical	Polymer	General Chemical Products
Headspace-GC/MS system	VOC in water	Flavor Water analysis	Residual solvents Water analysis	Evolved gas	Evolved gas
Purge & Trap GC/MS system	VOC in water	Water analysis	Water analysis		
SPME-GCMS system	Pesticides	Flavor Special componets		Evolved gas	Evolved gas
Pyrolizer-GC/MS system				Pyrolysis Evolved gas additives □	Evolved gas
Direct inlet system			Drug	Evolved gas	Evolved gas

Main Specifications

Specifications	GCMS-QP2010 Plus (Newly available)
GC	GC-2010
Ionization mode	EI/PCI/NCI
Ionization Voltage	10 ~200 eV
Emission Current	5uA ~250uA
Mass range	1.5 – 1090 m/z
Evacuation system	179L/s + 185L/s Dual inlet differential vacuum system
Maximum column flow rate	15mL/min
Sensitivity (OFN 1pg for EI) (BEN 100pg for CI) (OFN 100fg for NCI)	EI: S/N □ 160 CI: S/N □ 150 NCI: S/N □ 300
Scan Speed	10,000 u?sec
MS Interface Temperature	100 to 350°C
Ion source temperature	100 to 300°C
SIM	64 channels x 128 groups
AART	Standard
FASST	Standard

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