

# Theory and Utility of FT-NIR Spectroscopy


## An Introduction to Nicolet Antaris FT-NIR Analyzers

# MOLECULAR SPECTROSCOPY

**The study of the interaction of  
electromagnetic waves and matter**

(Includes UV, visible, NIR, IR, Raman etc.)

# Transitions Guide to the EM Spectrum



# SPECTRAL ABSORPTIONS

Radiowave  
Microwave  
IR

NMR  
Rotation of molecules, ESR  
Fundamental molecular  
vibrations

NIR

Overtones and combinations of  
Mid-IR

UV / Visible

Electronic transitions, energy of  
electrons raised to an excited  
state

X-Ray

Core electronic transitions in  
atom


# Vibrational Spectroscopic Regions

- Mid-IR Spectral Region
  - 4,000 - 400  $\text{cm}^{-1}$  (wavenumber)
  - 2,500 - 25,000 nm (wavelength of light)
- NIR Spectral Region
  - 12,820 - 4,000  $\text{cm}^{-1}$  (wavenumber)
  - 780 - 2,500 nm (wavelength of light)

To convert between nm and wavenumbers, simply take the inverse of the number in nm and multiply by  $10^7$ . For example, to convert 1000 nm into wavenumbers, take 1000, invert it and multiply by  $10^7$ . This gives  $10,000 \text{ cm}^{-1}$ .


# What Information Comes From the NIR ?

- Mid-IR spectroscopy makes use of *fundamental* vibrations
- Near-IR spectroscopy uses *Overtone* and *Combination Bands*
- Overtones and Combination Bands are types of vibrations related to the fundamental vibrations seen in the Mid-IR




# Nature of Overtones

- Musical instruments like guitars and pianos make use of fundamentals and overtones just like in spectroscopy
- Overtones occur at about 2 and 3 times the frequency of the fundamental vibration. Absorption intensity decreases with increasing overtones
- Band overlap increases with increasing overtones



# Fundamental and Overtones – A Quantum View





# Overtone

Occur Close to  
Integer Multiples Of Fundamental Bands.

For example: C-H Overtone Will Occur Near:

## First Overtone

$$2960\text{cm}^{-1} \text{ (C-H Stretch)} * 2 = 5920 \text{ cm}^{-1}$$

$$3378 \text{ nm (C-H Stretch)} * 1/2 = 1689 \text{ nm}$$


## Second Overtone

$$2960\text{cm}^{-1} \text{ (C-H Stretch)} * 3 = 8880 \text{ cm}^{-1}$$

$$3378 \text{ nm (C-H Stretch)} * 1/3 = 1126 \text{ nm}$$

# Combination Bands

- There are vibrations in the Near-IR other than overtones of fundamentals. These are called Combination Bands.
- COMBINATION Bands are the sum of several fundamentals from **different** vibrations. One photon of light excites two vibrations.
- Combination Bands are typically found at lower energies than overtones.




# Combination Bands

COMBINATION Bands Appear Near The Sum Of Two Or Three FUNDAMENTAL Bands

For Example: A C-H Combination Will Occur Near...

$$\begin{aligned} &2960\text{cm}^{-1} \text{ (C-H Stretch)} + 1460\text{cm}^{-1} \text{ (C-H Bend)} \\ &= 4420\text{cm}^{-1} \\ &= 2262\text{nm} \end{aligned}$$

# Example NIR Spectrum




# NIR Characteristics

- R-H groups have the strongest overtones due to hydrogen's degree of anharmonicity.
- O-H, N-H, C-H, S-H bonds etc. strong NIR absorbers
- H<sub>2</sub> and O<sub>2</sub> has no change in dipole moment, hence no NIR absorption
- R-H stretch or R-H stretch/bend form most NIR bands


# Typical Absorptions

— CH<sub>3</sub>  
 // CH<sub>2</sub>  
 \\\ CH  
 C=CH  
 Aromatic CH  
 // C=O  
 H<sub>2</sub>O  
 ROH  
 RNH<sub>2</sub>



# Typical FT-NIR Measurements

- Qualitative and quantitative analysis
  - *Identification*
  - *Qualification and quantitation*
  - *Determination of change*
  - *Process monitoring and control*
- Can measure all sample types



*Solids*

*Softgels*

*Films*

*Liquids*

*Grains*

*Suspensions*

*Tablets*

*Pastes*

*Syrups*

*Powders*

# Uses for NIR Spectroscopy

- NIR spectroscopy finds widespread use as an excellent QA/QC analysis tool in industrial applications:
  - *Pharmaceutical*
  - *Chemical*
  - *Petrochemical*
  - *Food and Beverage*
  - *Polymer*
  - *Agriculture*
- NIR is a proven and popular technique because it is rapid, accurate, non-destructive and can go through glass and plastic to analyze the material inside



# The FT-NIR Analysis Advantage

- Large pathlength
  - *Analyze more sample*
- Sampling *in situ* - through glass, packaging materials
  - *No sample preparation required*
- Accurate and precise
- Fast - Results in seconds
- No Hazardous reagents
- Non-destructive
- Remote sampling capabilities



# Regulatory Acceptance for Near IR



U.S. Pharmacopeia  
*The Standard of Quality<sup>sm</sup>*



- ASTM NIR Method E 1944
- USP Chapter <1119>
- European Pharmacopeia (Ph.Eur.) 2.2.40 1997  
(European Directorate for the Quality of Medicines)

# Types of Near IR Instrumentation


- Filter – Wavelength filters used to select specific wavelengths for analysis
- Dispersive – Gratings or prisms in combination with slits
- Fourier Transform (Nicolet Antaris) – Interferometer and internal HeNe laser

# Nicolet Antaris™ FT-NIR Analyzer

- A complete analyzer for the pharmaceutical, food & beverage, chemical, and polymer industries



# FT-NIR Spectrometer



# FT Technology for Near Infrared

- Based on the same interferometer technology used for mid-infrared and FT-Raman
- Utilizes all the advantages of FT based instruments
  - *Internal calibration is derived from He-Ne laser (precision =  $0.01\text{ cm}^{-1}$ )*
  - *High resolution ( $4\text{ cm}^{-1}$ )*
  - *All wavelengths are measured simultaneously*
  - *Higher energy throughput (larger apertures)*
  - *Excellent stability and reproducibility*
  - *Dynamic alignment*

# Limitation of Old Dispersive Technology


- Mechanically complicated
  - Grating filter and chopper
- Poor resolution
  - 16  $\text{cm}^{-1}$  or worse
  - Loss of spectral information
- Wavelength inaccuracy
  - Due to prism or grating
- Stray light effects
- Difficult to transfer methods due to wavelength inaccuracy



# Traditional NIR Problems

## Poor Wavelength Accuracy and Precision

- FT-NIR advantage: data repeatability and method reliability due to superior wavelength accuracy and precision
  - *Internal frequency reference with HeNe laser*






# Traditional NIR Problems

## Low Resolution



- Information can be hidden at lower resolution
  - *Need resolution to distinguish two closely related materials*
- FT-NIR advantage – simpler models with fewer calibration spectra possible with greater specificity



# Traditional NIR Problems

## Poor System Stability and Artifacts

- Stability – consistency in measurement results
- Artifacts – peaks or other anomalies in the spectrum not due to the sample
- Problems eliminated with Nicolet Antaris design



# Sampling with an FT-NIR Analyzer



Analyze • Detect • Measure • Control™

**Thermo**  
ELECTRON CORPORATION

# Measurement Techniques

- Transmission

- *Liquid analysis*
  - Quantitative
  - Qualitative
- *Films*
- *Tablets*
- *Semi-solids*



- Diffuse reflectance

- *Powders*
- *Solids*
- *Tablets*
- *Fibers*
- *Pastes*
- *Suspensions with high solids contents*
- *Opaque semi-solids*



# Nicolet Antaris FT-NIR Analyzer - Overview

**Integrating Sphere**

**Transmission**




**Fiber Optic**


Method Development Sampling System (MDS) shown above  
Dedicated configurations also available for specific analyses

# Nicolet Antaris FT-NIR Analyzers

- MDS System
  - All techniques on one system
- Other Dedicated Systems
  - RTS System
  - Liquid Analyzer
  - Solid Sampling System
  - Tablet Analyzer
  - Fiber Optic System
  - Fiber Optic Multiplexer



# Nicolet Antaris Optical Bench



- Rugged industrial platform
- Dynamic Alignment
- User-replaceable and pre-aligned consumables
- Modular design flexibility with stability of fixed beampaths
- Integrated validation (Built to meet or exceed USP monograph)

# Near IR Quartz Halogen Source

- External replacement with only 3 screws
- No wires
- Source change in seconds
- Pinned-in-place aluminum housing for lock-and-key fit
- High throughput






# Internal Referencing For All Collection Types

- Protected internal reference material simplifies background collection
  - *Spectralon*
  - *Diffuse Gold Flag*
- Eliminates errors due damaged reference materials
- Collect backgrounds with sample already loaded
- Easy to use and automated



# Reflection Sampling – Integrating Sphere

- Ideal collection environment for diffusely reflected light
- High efficiency collection
- Internal gold background material
- Sample simply needs to be placed on the instrument for data collection



# Liquid Sampling - Transmission Compartment

- Collect data on liquid, film or paper
- Accommodates vials, tubes, cuvettes and card holders
- Temperature control
- Automatic internal background position in transmission cell
- Reproducible positioning




# Fiber Optic Sampling – SabIR and 3<sup>rd</sup> Party Probes

- Transflectance or reflection sampling
- Fiber optic port uses standard SMA connectors
- Remote sampling capability for At-line, near-line, on-line
- Ideal for incoming raw material identification and in-process analysis




# Solid Transmission Sampling – Tablet Analyzer

- Thick, opaque tablets give quantitative results for
  - *Content uniformity*
  - *API or excipient concentrations*
  - *Physical properties like hardness*
- Integrating sphere reflection data gives tablet coating information
- Collect reflection and transmission data without removing the tablet
- Dedicated softgel detector




## Transmission



## Diffuse Reflection

# Antaris II

- Method Development System - All sampling types for NIR
  - *Transmission, tablets, softgels, diffuse reflection, fiber optics*
  - *Automatic sampling of solid dosage forms - 40 tablets at once*
  - *Simultaneous reflection and transmission measurement of tablets*
- Configurable base system with selection of sampling modules



# Remote Sampling – Fiber Optics – Antaris MX

- Simultaneous multi-channel measurement
- Simultaneous internal backgrounds
- No mechanical multiplexing
- Dual trigger probe capability
- Programmable analog/digital I/O options with RESULT Interface
- Spectral corrections adjust for fiber optic differences



Antaris MX




# Analyzer performance verification

- Traditionally done with probe pulled off process
- Internal reference channel can qualify operation while probe is inline
- Traceable reference standards for spectroscopic performance checks according to USP <1119>

## ParaLux™ illumination and DAQ system

All measurement points and reference data acquired at exactly the same time

Internal reference provides realtime referencing and internal qualification with traceable standard reference materials





# Fiber optic probes

- High quality probes with proprietary window mounts in custom configurations
- Specialty probes
  - *Fermentation – dual mode reflection and transmission*
  - *Retractable dryer probe with self cleaning tip and wash chamber*



# High Throughput Sampling – Autosampler RS and Multipro Autosampler

- Multiple samples with no operator interaction
- Eliminates sample-to-sample operator error
- Useful for analyzing tablets, softgels, solids or powders in vials
- Transfer already existing methods for seamless scale-up



# Heterogeneous Sampling – Sample Spinners

- Obtain spectra of heterogeneous samples like

*Polymer pellets*

*Grains*

*Coarse powders*

*Silage*

*Granular solids*

*Corn*

*Crystalline materials*

- Eliminate the need to collect multiple spectra on the same sample



# Antaris EX: A complete process analyzer

- Hazardous environments
- Integrated computing
- Integrated process communications
- Same sampling technology and advantages as MX
  - *Simultaneous measurements with internal referencing and qualification*
- Not a box that you have to put inside of another box!
- Complete new lineup of probes



# Hazardous environments

- Three EX models:
  - *EX-Z1 for Class I Division 1/Zone 1 Environments*
  - *EX-Z2 for Class I Division 2/Zone 2 Environments*
  - *EX-IP65 for*
- Enclosures are:
  - *Hose washable, Stainless steel, corrosion proof*
  - *Individually certified*
- Vortex cooling provides operating temperature range of -20 to 40 C
- Purge systems (different for each environment) linked to power of system





# On Blender Analysis



# Blending in Development – Multiple sizes

- Instrumentation must work with a wide range of blender sizes
- Issues of scale with equipment size and weight



# Monitoring the Blend with NIR


- By scanning the sample (pharmaceutical powder blend) through a window in the modified lid, the progression of mixing and concentration of components can be measured
- Data is triggered when bin is upside down (powder in contact with sampling window)
- Powered by rechargeable battery pack
- Diffuse reflection sampling
- Wireless communications to remote PC with RESULT software for control/DAQ
- Standard deviation measurement can be used to determine when mix is fully achieved



# The Spectroscopy Technology



Optical Bench



Optical Module Sub-Assembly




Spectrometer Assembly

- Source – Semiconductor based NIR tunable laser
  - Not a light bulb/globar – does not go “out” or fade*
  - Rated for a 25 year lifetime (Telcordia Spec.)*
- Wavelength selection - High resolution (2 nm / 8cm<sup>-1</sup>)
- Fabry-Perot tunable filter
  - Tuning range of hundreds of nm*
  - Scans full range in tens of milliseconds*
  - Insensitive to vibration*
- Internal Wavelength and Amplitude Reference
  - Ensures optimum wavelength calibration over time*



# Instrumentation and Sampling Interface

- **Detector: Single element InGaAs photodiode**
  - Optimized detector element sized to collect highest solid angle for maximum efficiency
- **Baseline s/n of 0.0002**
- **Wavelength ranges available**
  - 1350-1550 nm, 1550-1800 nm and 1800 – 1970 nm
- **Fast Scanning (~80 ms / scan)**
- **Sample sensing area of 40 mm (coincides with ~600 mg sample size)**
- **Smaller spot size available (~10 mm) for “intra-dosage” uniformity investigation**
- **Wireless communication (~15 meters)**
- **Battery operation (2 – 3 hours)**



# Summary

- On Blender Real-time analysis provides a rapid and calibration free method of getting blend information and performing experiments
- With the correct type of optical interface the instrumentation can be quickly moved between blender sizes and types for scale-up development
- Instrumentation is process 'friendly'
  - *Small*
  - *Lightweight*
  - *Easy to set up*
  - *Long battery life*
  - *Robust*



# RESULT Software Suite for the Nicolet Antaris Analyzer



Analyze • Detect • Measure • Control™

**Thermo**  
ELECTRON CORPORATION

# RESULT Software Suite for Near IR

- RESULT is the software package made for collecting Near IR data with the Nicolet Antaris Analyzer
- Provides
  - *Real time data acquisition*
  - *Sample prediction using developed models*
  - *Chemometrics tools for method development*
  - *Full suite of compliance tools*

# RESULT Provides Tools for all Users

- RESULT Operation allows users with no science training to run routine analyses
  - *Minimal training time and one-click use*
- RESULT Integration provides the tools to develop workflows with integrated SOPs for workroom deployment
  - *Design simple workflows like Raw Material ID*
  - *Create sophisticated experiments with logical gates, delays and repeats*

# RESULT Provides Tools for all Users

- TQ Analyst is a powerful chemometrics package
  - *Easy to use*
  - *Fully integrated with spectral collection software*
  - *No secondary programs needed*
- Administrator Mode
  - *Creates new users*
  - *Maintains system integrity*
  - *Sets workflow access privileges for all users*

# ValPro Instrument Qualification Platform


- Design qualification (DQ) documentation provided
- Detailed Installation Qualification (IQ) procedure
- USP-based Operational Qualification (OQ) tests provided
  - *Traceable & serialized internal validation wheel*
- Tools for Performance Qualification (PQ)





# 21 CFR Part 11 Compliance with RESULT

- Full Suite of Part 11 Tools
  - Full operation auditing
  - SOP enforcement and automation
  - Digital signatures
  - Security, access control, and privilege and workflow assignment



Key ID	Date	PC	Operator Name	Workflow Name	File Name	File Type	Item Name	Run Mode
275	08-24-2000 08:23:29	PC_DRAPER	Draper, Carla	Precision Constraint	C:\RESULT Data\Archive\precision 51731739-79C1-11D4-A31E- 0050DA5D91D9.htm	Report	statistics report	Production
273	08-24-2000 08:22:21	PC_DRAPER	Draper, Carla	Precision Constraint	C:\RESULT Data\Archive\ethanol\ethanol 51731736-79C1-11D4-A31E- 0050DA5D91D9.spa	Spectra	precision test	Production
272	08-24-2000 08:22:20	PC_DRAPER	Draper, Carla	Precision Constraint	C:\RESULT Data\Archive\ethanol\ethanol 51731735-79C1-11D4-A31E- 0050DA5D91D9.spa	Spectra	precision test	Production
271	08-24-2000 08:22:20	PC_DRAPER	Draper, Carla	Precision Constraint	C:\RESULT Data\Archive\ethanol\ethanol 51731734-79C1-11D4-A31E- 0050DA5D91D9.spa	Spectra	precision test	Production
270	08-24-2000 08:22:20	PC_DRAPER	Draper, Carla	Precision Constraint	C:\RESULT Data\Archive\ethanol\ethanol 51731733-79C1-11D4-A31E- 0050DA5D91D9.spa	Spectra	precision test	Production
269	08-24-2000 08:22:19	PC_DRAPER	Draper, Carla	Precision Constraint	C:\RESULT Data\Archive\ethanol\ethanol 51731732-79C1-11D4-A31E-	Spectra	precision test	Production

# Software Validation

- TQ Analyst Algorithms
  - *Fully Documented*
  - *Externally reproduced for verification in Excel*
  - *Demonstrates internal validation*

## Calibration

The following steps are performed in a PCR calibration for dimensions  $h = 1, 2, \dots, a$ :

1. Initialize a matrix  $\mathbf{U}$  that has  $r$  rows and has  $h = 0$  columns.
2. Increase  $h$  by 1 and select the column of  $\mathbf{X}$  (actually  $\mathbf{X}_h$  since it depends on the value of  $h$ ) with the greatest sum of squares. This is a first estimate for the principal component scores (or latent variables). Call this vector  $\mathbf{u}_h$ .  $\mathbf{u}_h$  is of size  $r \times 1$ .
3. Compute the squared norm of  $\mathbf{u}_h$ .

$$u_h^2 = \mathbf{u}_h' \mathbf{u}_h$$

4. Calculate the row vector  $\mathbf{b}_h$  as  $\mathbf{b}_h' = \mathbf{u}_h' \mathbf{X} / u_h^2$ .  $\mathbf{b}_h$  is of size  $p \times 1$ .

# Process Analytical Technology (PAT)

- PAT is a system for designing, analyzing, and controlling manufacturing through timely measurements (i.e., during processing) of critical quality and performance attributes of raw and in-process materials and processes with the goal of ensuring final product quality
- Seamless migration of FT-NIR technology from the laboratory to the manufacturing floor

# Time for a Paradigm Shift - PAT

- In-process testing that currently exists in pharmaceutical manufacturing consists of
  - *Physical measurements*
  - *Results from samples sent to the lab while processes were halted or continued “blindly”*
- Current testing is not real-time, nor “whole batch”
- Non-optimal conditions can still occur (over-blending for example)
- Testing and waiting for results is costly



# Typical FT-NIR PAT Applications



- Blending and mixing
- Reaction monitoring
- Fermentation
- Drying
- Post tableting tests




# NIR in the Manufacturing Process

- Traditionally, optimal mixing time was determined during process development and scale-up
- Avoids problems of incomplete or over-blending (de-blending)
- Improves efficiency of production processes and equipment usage
- Assures product uniformity
- Goals
  - *Establish adequacy of mixing*
  - *Eliminate thieving/process interference*
  - *Eliminate error and recipe approach*
  - *Accurate endpoint determination*



# Integration with manufacturing systems

- RESULT OPC (OLE for Process Control) Server facilitates powerful communication with manufacturing systems
- Send and receive commands automatically
- Trigger runs, endpoints, valves, etc.





# Communication is Essential to Process Control

- Communication system must be able to bridge
  - *Disparate technologies from various manufacturers*
  - *Compatibility with both high and low level interfaces*
  - *Opposing computer platforms*
  - *Geographic separation from target*
  - *Interface with business software to expand the process control loop*
- Data communication architecture must be inclusive, not exclusive
- Existing and modern facilities



# OPC in Process: Dryer Monitor



# Typical FT-NIR Pharmaceutical Applications

- Raw Material ID
- Moisture Content
- Content Uniformity
- Tablet API
- Coating Thickness
- Process Analytical Technology (PAT)



# Pharmaceutical Application: Raw material ID

- Apply NIR to inspect incoming raw materials quickly at the loading dock
- Equip the receiving area with an NIR to be used by technicians who are not trained as scientists
- Confirm the identity of each container of material that is received, without lab wet chemistry techniques that are time-consuming and expensive



# Library Raw Material Samples

- One lot each of the following substances was provided for the construction of the raw material library

d-Glucose  
d-Fructose  
Sucrose  
d-Mannitol  
d-Sorbitol  
 $\alpha$ -d-Lactose Monohydrate  
Acetylsalicylic Acid  
Acetaminophen  
L-Ascorbic Acid  
Citric Acid



# Validation Samples

d-Glucose - Second lot from a different manufacturer

d-Fructose - Second lot from a different manufacturer

Acetaminophen - Second lot from a different manufacturer

**Salicylic Acid** - Degradation product of Acetylsalicylic Acid

**2-acetamidophenol** - Constitutional isomer of acetaminophen

**$\alpha$ -d-Lactose anhydrous** - Dehydrated form of library material

Black items are positive challenge samples (Should Pass)

**Red** items are negative challenge samples (Should Fail)

## External Challenge Validation Results

<u>Compound</u>	<u>Score</u>	<u>Result</u>
Glucose (Lot #2)	99.9	Pass
Fructose (Lot #2)	98.0	Pass
Acetaminophen (Lot #2)	100.0	Pass
Salicylic Acid	55.2	Fail
2-Acetamidophenol	19.3	Fail
Lactose Anhydrous	68.9	Fail

## FT-NIR Success stories

- Customer had a plant with diverse array of incoming raw materials that needed to be classified accurately.
- Implemented NIR testing for over 200 materials
- Based upon technician time savings alone, paid back the investment in about 7 months
- Needed to identify a small number of materials at very high volume composing the bulk of their test load
- Implemented 100 % inspection for ID and qualification
- Dropped their testing time on each incoming shipment of raw material by 68 %



# Lyophilized Material Analysis

- Lyophilized materials are an excellent application for analysis by NIR
- Packaging cannot be opened without corrupting the product inside
- Can apply to food, pharmaceutical, fine chemical or biotech industries






# Current Testing for Lyophilized Products

- Current testing methods like Karl Fisher titration, LOD or GC
  - *cumbersome, time consuming and costly*
  - *destructive*
- NIR does not need solvents or carrier gases
- NIR calibrations for lyophilized materials have been shown to be as accurate as the reference method
- Other components can also be measured using the same NIR spectrum (Buffer, API, Cryoprotectant)


# Case Study For Lyophilized Sample Analysis

- Primary component of interest in any lyophilized materials is moisture (0.5% - 2.0% typically)
- Sample is Thrombin, a topical coagulant used in the medical and dental fields
- Thrombin loses its efficacy if exposed to high heat or wet conditions. Current application will look at moisture and potency in Thrombin


# Moisture First Overtone Region




# Close-up of First Moisture Overtone Region



# Calibration Curve for Moisture



# Calibration for Potency



# Quantification of Active Pharmaceutical Ingredient (API) in Tablets


## Show Concentration Differences in Tablets

- Distinguish different clinical tablet formulations by amount of active ingredient
- Finished product is a tablet and manufacturing protocols need independent verification of amount of active ingredient
- Understand the differences between transmission and reflection analysis of tablets







# Reflectance Spectra for Clinical Tablets




# Tablet Reflection Results



# Tablet Transmission Results



# API Quantification Can Also Be Done Using Qualitative Methodology - Discriminant Analysis



- Principal component scores plot reveals quantitative method may not be necessary
- The required information is obtained with excellent results using qualitative analysis (Discriminant)

# NIR in Pharmaceuticals - Summary

- NIR has proven useful in the following applications:
  - *Raw Material ID*
  - *Raw Materials - Particle Size*
  - *Analysis of Parenteral Drugs - Lyophilized Materials*
  - *Quantification of API in Tablets*
  - *Blend Homogeneity*
  - *Lactose Polymorph Screening*
  - *SoftGel Analysis*
  - *Enzyme Activity*

# Near IR in the Chemical Industry

- Near IR analysis has proven worth across the chemical industry
  - *Fine chemicals*
  - *Polymers*
  - *Petrochemicals*
  - *Paints, coatings and adhesives*
  - *Soaps and cleaners*
  - *Synthetic materials*
- At-line, near-line and on-line



# Polymerization Monitoring: Acid Value and Hydroxide

Analyze • Detect • Measure • Control™

**Thermo**  
ELECTRON CORPORATION


## HoVal Analysis in NIR - Purpose

- HoVal is short for **H**ydroxyl **V**alue, an important metric in polymerization reactions
- Measures both starting materials (Acid and OH) simultaneously from the same spectra
- Condensation polymers (Polyesters) show decrease in acid and alcohol concentration with time
- Use NIR to monitor this decrease will show when polymerization is finished




# Experimental


- Samples were placed into 6 mm disposable vials
- Analysis using Nicolet Antaris Transmission module



# Spectral Data for HoVal Analysis




# NIR Replaces Two Wet Chemical Methods



Hydroxyl Value Correlation

Acid Number Correlation



# Plasticizer Content in PVC

# Plasticizer Content in PVC

- Plasticizers are common additives to polymers.
- The quantity of plasticizer determines the properties of the polymer.
- The plasticizer contents of polymers differ depending on the application.
- Primary methods often involve tedious and time-consuming extractions.

# Samples

- Samples of polyvinyl chloride (PVC) with varying amounts of dioctylphthalate (DOP) as plasticizer
- Three types of samples:
  - *Opaque PVC plates with 5 - 50% DOP*
  - *Transparent PVC plates with 5 - 50% DOP*
  - *PVC Films with 9 - 40% DOP*




# FT-NIR Measurements

- *As is* measurements (no sample preparation or modification).
- Data collected from 4,000  $\text{cm}^{-1}$  to 10,000  $\text{cm}^{-1}$
- 4  $\text{cm}^{-1}$  resolution.
- 90 scans
- One (1) minute collection time.
- All measurements done on Integrating Sphere



# Non-Transparent Plates






# Calibration Parameters

	Factors	Data Format	Smoothing	Region (cm <sup>-1</sup> )	Pathlength
Opaque	4	2 <sup>nd</sup> Derivative	Norris 9,0	8982-3920	MSC
Transparent	4	2 <sup>nd</sup> Derivative	Norris 5,5	4772-4520 5935-5322	MSC
Films	3	2 <sup>nd</sup> Derivative	Norris 5,5	4125-3946 4296-4263	MSC

# Non-Transparent Plates



# Calibration Results

	Calibration		Cross Validation	
	<u>Corr. Coeff.</u>	<u>RMSEC</u>	<u>Corr. Coeff.</u>	<u>RMSECV</u>
Opaque	0.999	0.375	0.995	1.49
Transparent Plates	0.999	0.591	0.996	1.45
Films	0.999	0.174	0.999	0.407

# Conclusions on NIR of Plasticizer Content

- Additives such as plasticizers can be measured in various forms of plastic materials
- These materials can be finished or unfinished
  - *Measurements can be in-process or QC measurements*
- The analyses can be done easily, rapidly and non-destructively
- Current methodology for extraction of matrix elements in polymers is extremely time consuming

# Near IR and Food

- Near IR analysis of food and food products has found extensive use in all segments of the industry
  - *Processed foods*
  - *Beverages*
  - *Grain*
  - *Fruit*
  - *Meat*
- Common components like protein, moisture, ash, starch and fat are measured simultaneously without destroying the sample



# FT-NIR As A Process Analytic in Food

- Develop methods for any type of food product
- High resolution advantage of FT technology keeps number of standards low
- Many different sampling accessories for food products
  - *Sample spinner*
  - *SablR Fiber Optic Probe*
  - *Cheese and Viscous Liquid Samplers*
  - *Softgel detector (Gum)*
- OPC Output for analyzer incorporation into Process Control Systems

# NIR Analysis in the Food Industry

- Moisture and Protein in Barley
- Viscous Liquid Analysis
- Oil in Rice Powder
- Amylose Levels in Corn
- Multi-component Analysis of Wine
- Fat and Dry Matter in Cheese
- Discrimination of Jellybeans
- Choline Chloride Levels in Silage

# Nicolet Antaris Product Line – Complete Solution

- FT-NIR Nicolet Antaris Spectrometer Platform



- RESULT Analyzer Software Platform




- TQ Analyst Chemometric Software Platform



- ValPro Validation/Qualification Platform



- Support “Pyramid of Support” products and services





# The FT-NIR Advantage

The Bottom Line:  
**Save Time and Money.**

**Use FT-NIR to its maximum capability.**